



CX41-32C02详细配置清单

E/C Code	Description	quantity
CX41-32C02 显微镜镜体	Biological Microscope (Made in Philippines)	1
CX41RF 三目镜筒	Microscope frame with quintuplerevolving nosepiece, mechanicalstage (right handle), 30W illuminator, filter 32.5C-2. Corresponded to WEEE regulation	1
U-CTR30-2 附加透镜	Trinocular tube, FN20, fixed light pass, Bi/Photo: 50/50	1
CX-AL 明场聚光镜	Attachment lens for condenser CH3-CD, CH3-CDP, CX-SLC	1
CH3-CD 6V30W 灯泡	Brightfield condenser, N.A. 1.25, with A.S. vernier scale	1
6V30WHAL 电源线	30W Halogen bulb (Philips No. 5761)	2
UYCP 防尘罩	Power cord	1
COVER-015 平场消色差物镜	Dust cover Type 015 for CX21/31/41 (L 450mm x H 460mm)	1
PLCN4X	Plan achromat objective 4X/0.1, WD 18.5	1
PLCN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLCN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
PLCN100XO 目镜	Plan achromat objective 100X/1.25, WD 0.13 (spring, oil)	1
WHB10X 镜油	Wide field eyepiece 10X, FN20	2
OIL-8CC	Immersion oil 8cc	1