



BX53T-32P02 详细配置清单

E/C Code	Description	Q
BX53T-32P02	Research Microscope	1
BX53机架	Microscope frame for transmitted and reflected microscopy with built-in filter holders for four filters (LBD, ND6, ND25 and GF) , including hexagonal driver ABxxxx and immersion oil 8cc Corresponded to WEEE and RoHS regulations	1
BX53F	Microscope frame for transmitted and reflected microscopy with built-in filter holders for four filters (LBD, ND6, ND25 and GF) , including hexagonal driver ABxxxx and immersion oil 8cc Corresponded to WEEE and RoHS regulations	1
三目观察筒	Trinocular tube	1
U-TR30-2	Trinocular tube	1
D 型6孔物镜转盘	Sextuple revolving nosepiece with a slot for analyzer or DICslider	1
U-D6RE	Sextuple revolving nosepiece with a slot for analyzer or DICslider	1
右手低位载物台	Ceramic surface mechanical stage with right-hand low drivecontrol (long type)	1
U-SVRB-4	Ceramic surface mechanical stage with right-hand low drivecontrol (long type)	1
左手双标厚型标本夹	Specimen holder for BX stage, left hand, for two specimens, thick type	1
U-HLDT-4	Specimen holder for BX stage, left hand, for two specimens, thick type	1
摇入摇出型聚光镜	Swing-out condenser N.A. 0.9-0.16	1
U-SC3	Swing-out condenser N.A. 0.9-0.16	1
100W 卤素灯灯室	Lamp house for 100W halogen	1
U-LH100-3	Lamp house for 100W halogen	1
12伏100瓦长寿命卤素灯泡	Halogen bulb (Philips No. 7724)	2
JC12V100WHAL-L	Halogen bulb (Philips No. 7724)	2
电源线	Power cord	1
UYCP	Power cord	1
BX3镜臂	Arm for BX53F and BX63F frame	1
BX3-ARM	Arm for BX53F and BX63F frame	1
防尘罩018	Dust cover Type 018 for BX2 (L 510mm x W 320mm x H 530 mm)	1
COVER-018	Dust cover Type 018 for BX2 (L 510mm x W 320mm x H 530 mm)	1
平场消色差物镜	Plan achromat objective 4X/0.1, WD 18.5	1
PLN4X	Plan achromat objective 4X/0.1, WD 18.5	1
PLN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
PLN100XO	Plan achromat objective 100X/1.25, WD 0.15 (spring, oil)	1
10倍宽视场目镜	Widefield eyepiece 10X	1
WHN10X	Widefield eyepiece 10X	1

10倍宽视场可调目镜

WHN10X-H

Widefield eyepiece 10X, focusable

1
