



**BX43-32P01HAL 详细配置清单**

<b>E/C Code</b>	<b>Description</b>	<b>quantity</b>
<b>BX43-32P01HAL</b>	<b>Laboratory Microscope</b>	<b>1</b>
BX43机架		
	Microscope frame for transmitted microscopy, including hexagonal driver ABxxxx, immersion oil 8cc Corresponded to RoHS, WEEE regulation	1
BX43F		
三目观察筒		
U-TR30-2	Trinocular tube	1
一般型5孔转盘		
U-5RE-2	Quintuple Revolving Nosepiece	1
右手低位载物台		
U-SVRB-4	Ceramic surface mechanical stage with right-hand low drivecontrol (long type)	1
30瓦卤素灯适配器		
U-LS30ADP	Adapter for attaching LS30 to BX43/43C frame, with connection cable Corresponded to RoHS, WEEE regulations	1
左手双标厚型标本夹		
U-HLDT-4	Specimen holder for BX stage, left hand, for two specimens, thick type	1
摇入摇出型聚光镜		
U-SC3	Swing-out condenser N.A. 0.9-0.16	1
30瓦卤素灯电源		
	Power supply for 30W halogen lamp with output cable (0.64m) (Made in Philippines) Corresponded to WEEE regulation	1
TL-4		
30瓦卤素灯灯座		
U-LS30-5	Lamp socket	1
6伏30瓦卤素灯		
6V30WHAL	30W Halogen bulb (Philips No. 5761)	2
日光平衡滤色片		
45-LBD-IF	Interference light balance daylight filter, 45mm dia.	1
电源线		
UYCP	Power cord	2
防尘罩018		
COVER-018	Dust cover Type 018 for BX2 ( L 510mm x W 320mm x H 530 mm )	1
平场消色差物镜		

PLN4X	Plan achromat objective 4X/0.1, WD 18.5	1
PLN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLN20X	Plan achromat objective 20X/0.4, WD 1.2 (spring)	1
PLN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
PLN100XO	Plan achromat objective 100X/1.25, WD 0.15 (spring, oil)	1
10倍宽视场目镜		
WHN10X	Widefield eyepiece 10X	1
10倍宽视场可调目镜		
WHN10X-H	Widefield eyepiece 10X	1
8毫升镜油		
OIL-8CC	Immersion oil 8cc	1

---