



CX31-32C02详细配置清单

E/C Code	Description	quantity
CX31-32C02 显微镜镜体	Biological Mirroscope (Made in Philippines)	1
CX31RTSF	Microscope frame with trinocular tube, a pair of eyepiece 10X (FN20), quadruple revolving nosepiece, mechanicalstage (right handle), abbe condenser, Plan objective 4X, 10X, 40X, AS, blue filter, 30W illuminator, 6V30W halogen bulb. Corresponded to WEEE regulation	1
6V30W 灯泡		
6V30WHAL	30W Halogen bulb (Philips No. 5761)	1
电源线		
UYCP	Power cord	1
防尘罩		
COVER-015	Dust cover Type 015 for CX21/31/41 (L 450mm x H 460mm)	1
平场消色差物镜		
PLCN100XO	Plan achromat objective 100X/1.25, WD 0.13 (spring, oil)	1
镜油		
OIL-8CC	Immersion oil 8cc	1