



CX41-12C02详细配置清单

E/C Code	Description	quantity
CX41-12C02	Biological Microscope (Made in Philippines)	1
显微镜镜体		
CX41RF	Microscope frame with quintuplerevolving nosepiece, mechanicalstage (right handle), 30W illuminator, filter 32.5C-2. Corresponded to WEEE regulation	1
双目镜筒		
U-CBI30-2	Binocular tube, FN20	1
附加透镜		
CX-AL	Attachment lens for condenser CH3-CD, CH3-CDP, CX-SLC	1
明场聚光镜		
CH3-CD	Brightfield condenser, N.A. 1.25, with A.S. vernier scale	1
6V30W 灯泡		
6V30WHAL	30W Halogen bulb (Philips No. 5761)	2
电源线		
UYCP	Power cord	1
防尘罩		
COVER-015	Dust cover Type 015 for CX21/31/41 (L 450mm x H 460mm)	1
平场消色差物镜		
PLCN4X	Plan achromat objective 4X/0.1, WD 18.5	1
PLCN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLCN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
PLCN100XO	Plan achromat objective 100X/1.25, WD 0.13 (spring, oil)	1
目镜		
WHB10X	Wide field eyepiece 10X, FN20	2
镜油		
OIL-8CC	Immersion oil 8cc	1