



BX43-12P01 详细配置清单

E/C Code	Description	Q
BX43-12P01	Laboratory Microscope	1
BX43机架		
	Microscope frame for transmitted microscopy, including hexagonal driver ABxxxx, immersion oil 8cc Corresponded to RoHS, WEEE regulation	1
BX43F		
通用双目镜筒		
U-BI30-2	Binocular tube	1
BX43编码型5孔物镜转盘		
BX43-5RES	Coded quintuple nosepiece for BX43	1
右手低位载物台		
U-SVRB-4	Ceramic surface mechanical stage with right-hand low drivecontrol (long type)	1
左手双标厚型标本夹		
U-HLDT-4	Specimen holder for BX stage, left hand, for two specimens, thick type	1
万能低倍聚光镜		
U-LC	Low magnification condenser	1
LED 灯室		
U-LHLEDC	Lamp house for LED with connecting cable Corresponded to RoHS, WEEE regulations	1
电源线		
UYCP	Power cord	1
防尘罩018		
COVER-018	Dust cover Type 018 for BX2 (L 510mm x W 320mm x H 530 mm)	1
平场消色差物镜		
PLN2X	Plan achromat objective 2X/0.06, WD 5.8	1
PLN4X	Plan achromat objective 4X/0.1, WD 18.5	1
PLN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLN20X	Plan achromat objective 20X/0.4, WD 1.2 (spring)	1
PLN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
10倍宽视场目镜		
WHN10X	Widefield eyepiece 10X	1
10倍宽视场可调目镜		
WHN10X-H	Widefield eyepiece 10X, focusable	1