

**BX46-72P05详细配置清单**

E/C Code	Description	quantity
BX46-72P05	Laboratory Microscope	1
BX46机架		
	Microscope frame for transmitted microscopy, including	
BX46F	hexagonal driver ABxxxx, immersion oil 8cc Corresponded to	1
	RoHS, WEEE regulation	
通用可倾斜双目观察筒		
U-TBI3	Tilting binocular tube	1
右手低位载物台		
U-SVRB-4	Ceramic surface mechanical stage with right-hand low	1
	drivecontrol (long type)	
左手单玻片夹片器		
U-HLS-4	Specimen holder for BX stage, left hand, for single specimen	1
橡胶帽		
U-SHG	Rubber caps for BX stage knobs, including for X-axis knob (1	1
	pc.) and for Y-axis knob (1 pc.) , thin type	
LED 灯室		
U-LHLEDC	Lamp house for LED with connecting cable Corresponded to	1
	RoHS, WEEE regulations	
电源线		
UYCP	Power cord	1
防尘罩018		
COVER-018	Dust cover Type 018 for BX2 (L 510mm x W 320mm x H 530	1
	mm)	
平场消色差物镜		
PLN4X	Plan achromat objective 4X/0.1, WD 18.5	1
PLN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
10倍宽视场目镜		
WHN10X	Widefield eyepiece 10X	1
10倍宽视场可调目镜		
WHN10X-H	Widefield eyepiece 10X, focusable	1