



### CKX41-A22PHP 详细配置清单

Units	Description	Qty
显微镜镜体	CKX41 frame, including 30W lampsocket AR3553, insert plate AQ1180, allen wrench GA3720, tag AD0764, 6V30WHAL (2pcs) and AB4X10SA screw Corresponded to WEEE regulation	1
CKX41SF		
可倾斜式双目镜筒	Tilting binocular tube for CKX41, including eyepieces AR6118 (2 pcs.)	1
CKX-TBI		
载物台延长板	Stage extension plate	1
CK2-SS		
机械式载物台	Attachable mechanical stage with right-hand vertical low drive controls (including Petri-dish holder, Terasaki- plate holder, and Slide glass holder)	1
CK40-MVR		
玻璃插板	Glass stage insert plate for IX50/70, CK40F and CKX41SF	1
CK40-CPG30		
相差插板	Phase contrast slider, pre- centered	1
IX2-SLP		
电源线	Power cord	1
UYCP		
防尘罩	Dust cover Type 012 for CKX, CK (L 300mm x W 440mm x H 440mm)	1
COVER-012		
万能平场半复消色差相差物镜	U plan semi apochromat phase (pre-center) objective 4X/0.13, W.D. 17.0	1
UPLFLN4XPHP		
消色差相差物镜	Acromat phase (pre-center) objective 10X/0.25, W.D. 8.8mm	1
CACHN10XPHP		
长工作距离消色差相差物镜	Long working distance achromat phase contrast (pre-center) objective 20X/0.4, WD 3.2	1
LCACHN20XPHP		
LCACHN40XPHP	Long working distance achromat phase contrast (pre-center) objective 40X/0.55, WD 2.2	1
日光平衡滤色片	Interference light balance daylight filter, 45mm dia.	1
45-LBD-IF		