



CKX31-A12PHP 详细配置清单

E/C Code	Description	quantity
CKX31-A12PHP 显微镜镜体	Inverted Microscope (Made in Philippines)	1
CKX31SF 载物台延长板	CKX31 frame with binocular tube and eyepieces (2pcs) , including 30W lamp socket AR3553, 35mm dia. Petri dish holder AD0675, allen wrench GA3720, tag AD0764, 6V30W HAL (2pcs) and AB4X10SA screw Corresponded to WEEE regulation	1
CK2-SS 机械式载物台	Stage extension plate	1
CK40-MVR 相差插板	Attachable mechanical stage with right-hand vertical low drive controls (including Petri-dish holder, Terasaki- plate holder, and Slide glass holder)	1
IX2-SLP 电源线	Phase contrast slider, pre- centered	1
UYCP 防尘罩	Power cord	1
COVER-012 万能平场半复消色差 相差物镜	Dust cover Type 012 for CKX, CK (L 300mm x W 440mm x H 440mm)	1
UPLFLN4XPHP 消色差相差物镜	U plan semi apochromat phase (pre-center) objective 4X/0.13, W.D. 17.0	1
CACHN10XPHP 长工作距离消色差相 差物镜	Acromat phase (pre-center) objective 10X/0.25, W.D. 8.8mm	1
LCACHN20XPHP 日光平衡滤色片	Long working distance achromat phase contrast (pre-center) objective 20X/0.4, WD 3.2	1
45-LBD-IF	Interference light balance daylight filter, 45mm dia.	1