



IX71-F22FL/PH 详细配置清单

Units	Description	Qty
IX71-F22FL/PH 显微镜镜体	Inverted Microscope	1
IX71S1F3	Microscope stand with Light ON/OFF switch, sextuple revolving nosepiece, magnification changer (1X and 1.6X) and 2-step light path selector (observation:left sideport=100:0, 0:100), corresponded to near-infrared(700-1000nm), relay lens, CA1.6X, prism, dust proof glass-AW multi coated	1
中性灰度滤色片		
45-ND6	Neutral density filter, 45mm dia.	1
45-ND25	Neutral density filter, 45mm dia.	1
32ND25	Neutral density filter	1
32ND6	Neutral density filter	1
双目镜筒		
U-BI90CT 玻璃插板	Binocular tube with CT	1
IX2-GCP	Glass insert plate with 25mm dia. Opening, including rubber bands for objective lens (5 pcs; Red (AD9322), Yellow (AD9323), Green (AD9324), Blue (AD9325), White (AD9326), Light Blue (AD9327))	1
机械式载物台		
IX2-SFR	2-way mechanical stage with flexible right-hand vertical low drive controls, including stage insert plate AQ1180 and screws AB5X12SA (4 pcs)	1
超长工作距离聚光镜		
IX-ULWCD	Ultra long working distance condenser for phase contrast (N.A.0.30/W.D.73mm), four position (three positions for optical element and one for bright field), including IX-PHLU and IX-PH2U	1
相差环板		
IX-PHCU 灯柱	Phase contrast ring slit for IX-ULWCD	1
IX2-ILL100	Illuminator pillar (tiltable) for 100W HAL with condenser holder and 4 filter holders (45mm dia.), filed stop, including diffusion filter LP4539, rubber cap AD4797, cord clamp DO1492 and fixing screws AB5X35SA (4 pcs)	1
100W 汞灯灯箱		
U-LH100HG	Lamp house for 100W mercury burner with collector lens and connecting cable	1

100W 卤素灯灯箱		
U-LH100L-3	Lamp house for 100W halogen including long connection cord Corresponded to RoHS, WEEE regulations	1
聚焦透镜调节手柄		
U-CLA	Collector lens adjustment knob for U-ULH	1
供电箱		
TH4	Power supply	1
汞灯供电箱		
U-RFL-T	Power supply unit for Mercury lamp meeting the RoHS regulation	1
100W 卤素灯灯泡		
JC12V100WHAL-L	Halogen bulb (Philips No. 7724)	2
电源线		
UYCP	Power cord	2
荧光激发块		
U-MWB2	Mirror Unit for fluorescence microscopy, wide band Blue excitation, exciter filter BP450-480, dichroic beamsplitter DM500, barrier filter BA515	1
U-MWG2	Mirror Unit for fluorescence microscopy, wide band Green excitation, exciter filter BP510-550, dichroic beamsplitter DM570, barrier filter BA590	1
U-MWU2	Mirror Unit for fluorescence microscopy, wide band UV excitation, exciter filter BP330-385, dichroic beamsplitter DM400, barrier filter BA420	1
对中屏		
U-CST	Centering specimen with target for fluorescence microscopy	1
荧光激发块转盘		
IX2-RFAC	Fluorescence mirror unit cassette (6 position) with shutter and shockless mechanism, including color beads for indication (Blue: AD9555 (6 pcs), Red: AD9556 (6 pcs) and Green: AD9557 (6 pcs))	1
荧光照明装置		
IX2-RFA	Fluorescence illuminator (340nm applicable) with filter slider (32mm dia.), including UV light protection plate AW0243	1
防尘罩		
COVER-020	Dust cover Type 020 for IX2, MVX10 (L 620mm x W 250mm x H 660 mm)	1
100W 汞灯		
USH103OL	100W mercury burner for U-LH100HGAPO, U-LH100HG	2
平场消色差相差物镜		
CPLN10XPH	Plan achromat phase objective 10X, NA 0.25, WD 10	1
长工作距离消色差相 差物镜		
LCACHN20XPH	Long working distance achromat phase objective 20X/0.4, WD	1

	3.2	
LUCPLFLN40XPH	Long working distance plan semiapochromat phase objective 40X, NA 0.6, WD 3.0-4.2 (c.c.0-2)	1
目镜		
WHN10X	Widefield eyepiece 10X	1
WHN10X-H	Widefield eyepiece 10X, focusable	1
绿色滤光片		
43IF550-W45	Interference green contrast filter 45mm dia.	1
