



IX51-A12PH 详细配置清单		
E/C Code	Description	quantity
IX51-A12PH	Inverted Microscope	1
显微镜镜体		
IX51S8F3	Microscope stand with sextuple revolving nosepiece and two step light path selector (observation:left side port= 100:0, 20:80), corresponded to near-infrared(700-1000nm), relay lens, prism, dust proof glass-AW multi coated	1
双目镜筒		
U-BI90	Binocular tube	1
机械式载物台		
IX2-SFR	2-way mechanical stage with flexible right-hand vertical low drive controls, including stage insert plate AQ1180 and screws AB5X12SA (4 pcs)	1
相差插槽		
IX2-SL	Phase contrast slider, centrable, with phase ring for 4X, phase ring for 10X/20X and hole for phase ring for 40X (IX2-SLPH2) or filter (45mm dia.), including centering tools AQ2302 (2 pcs)	1
相差环板		
IX2-SLPH2	Phase contrast ring slit	1
IX2-SLPHC	Phase contrast ring slit	1
灯柱		
IX2-ILL30	Illumination pillar for 30W HAL with one filter holder (45mm dia.), including fixing screws AB5X35SA (4 pcs)	1
卤素灯供电箱		
TL-4	Power supply for 30W halogen lamp with output cable (0.64m) Corresponded to WEEE regulation	1
30W 卤素灯插座		
U-LS30-3-5	Lamp socket	1
30W 卤素灯灯泡		
6V30W HAL	30W Halogen bulb (Philips No. 5761)	2
电源线		
UYCP	Power cord	1
防尘罩		
COVER-020	Dust cover Type 020 for IX2, MVX10 (L 620mm x W 250mm x H 660 mm)	1
万能平场半复消色差相差物镜		
UPLFLN4XPH	U plan semi apochromat phase objective 4X/0.13, WD 17	1
平场消色差相差物镜		
CPLN10XPH	Plan achromat phase objective 10X, NA 0.25, WD 10	1

长工作距离消色差相差物镜			
LCACHN20XPH	Long working distance achromat phase objective 20X/0.4, WD 3.2		1
LCACHN40XPH	Long working distance achromat phase objective 40X/0.55, WD 2.2		1
目镜			
WHN10X	Widefield eyepiece 10X		1
WHN10X-H	Widefield eyepiece 10X, focusable		1
对中望远镜			
U-CT30	Centering telescope		1
绿色滤光片			
43IF550-W45	Interference green contrast filter 45mm dia.		1