



CKX41-F32FL 倒置荧光显微镜详细配置清单

Units	Description	Qty
CKX41-F32FL 显微镜镜体	Inverted Microscope (Made in Philippines)	1
CKX41SF 中性灰度滤色片	CKX41 frame, including 30W lampsocket AR3553, insert plate AQ1180, allen wrench GA3720, tag AD0764, 6V30WHAL (2pcs) and AB4X10SA screw Corresponded to WEEE regulation	1
32ND25 三目镜筒	Neutral density filter	1
U-CTR30-2 载物台延长板	Trinocular tube, FN20, fixed light pass, Bi/Photo: 50/50	1
CK2-SS 机械式载物台	Stage extension plate	1
CK40-MVR 玻璃插板	Attachable mechanical stage with right-hand vertical low drive controls (including Petri-dish holder, Terasaki- plate holder, and Slide glass holder)	1
CK40-CPG30 汞灯灯室	Glass stage insert plate for IX50/70, CK40F and CKX41SF	1
U-LH50HG 汞灯供电箱	Lamp housing for 50W mercury	1
U-RFLT50 电源线	Power supply for mercury	1
UYCP 荧光照明装置	Power cord	2
CKX-RFA*	Reflected light fluorescence illuminator equipped with field stop, 3-position fluorescence slider (with B excitation filter set and G excitation filter set) and 3-position (ND, empty and shutter) ND slider, including UV shield plate AV7728, operation knob AE1678, allen wrenches AC7359 (2 pcs) and AB4X10SA screws (3 pcs)	1
防尘罩		
COVER-012	Dust cover Type 012 for CKX, CK (L 300mm x W 440mm x H 440mm)	1
平场消色差物镜		
PLCN4X	Plan achromat objective 4X/0.1, WD 18.5	1
PLCN10X 长工作距离万能平场	Plan achromat objective 10X/0.25, WD 10.6	1

半复消色差物镜		
LUCPLFLN40X	Long working distance Plan SemiApochromat objective 40X, NA 0.6, WD 2.7-4.0 (c.c.0-2)	1
宽视野目镜		
WHB10X	Wide field eyepiece 10X, FN20	2
汞灯灯泡		
HBO50W-AC	50W Mercury burner for U-LH50HG	2
日光平衡滤色片		
45-LBD-IF	Interference light balance daylight filter, 45mm dia.	1
中性灰度滤色片		
43ND-25W45	Neutral density filter, 45mm dia.	1
