



**IX51-F32PH 详细配置清单**

<b>E/C Code</b>	<b>Description</b>	<b>quantity</b>
IX51-F32PH 显微镜镜体	Inverted Microscope	1
IX51S8F3 日光平衡滤光片	Microscope stand with sextuple revolving nosepiece and two step light path selector (observation:left side port= 100:0, 20:80), corresponded to near-infrared(700-1000nm), relay lens, prism, dust proof glass-AW multi coated	1
45-LBD-IF 中性灰度滤色片	Interference light balance daylight filter, 45mm dia.	1
45-ND6	Neutral density filter, 45mm dia.	1
45-ND25	Neutral density filter, 45mm dia.	1
IX-SPT 单口镜筒	Single port tube for U-TR30H trinocular tube, and IX70/50 side port	1
IX-ATU 中间镜筒	Intermediate tube attachment for 45 degree light path reflection	1
U-TR30H-2 三目镜筒	Wide field trinocular tube, ecofriendly material is used for the relay lens	1
IX2-SFR 机械式载物台	2-way mechanical stage with flexible right-hand vertical low drive controls, including stage insert plate AQ1180 and screws AB5X12SA (4 pcs)	1
IX-ULWCD 超长工作距离聚光镜	Ultra long working distance condenser for phase contrast (N.A.0.30/W.D.73mm), four position (three positions for optical element and one for bright field), including IX-PHLU and IX-PH2U	1
IX-PH1U 相差环板	Phase contrast ring slit for IX-ULWCD	1
IX2-ILL100 灯柱	Illuminator pillar (tiltable) for 100W HAL with condenser holder and 4 filter holders (45mm dia.), filed stop, including diffusion filter LP4539, rubber cap AD4797, cord clamp DO1492 and fixing screws AB5X35SA (4 pcs)	1
U-LH100L-3 100W 卤素灯灯箱	Lamp house for 100W halogen including long connection cord Corresponded to RoHS, WEEE regulations	1
TH4-HS 供电箱手动面板	Hand switch for TH4	1

供电箱			
TH4	Power supply		1
100W 卤素灯灯泡			
JC12V100WHAL-L	Halogen bulb (Philips No. 7724)		2
电源线			
UYCP	Power cord		1
防尘罩			
COVER-020	Dust cover Type 020 for IX2, MVX10 (L 620mm x W 250mm x H 660 mm)		1
万能平场半复消色差 相差物镜			
UPLFLN4XPH	U plan semi apochromat phase objective 4X/0.13, WD 17		1
UPLFLN10X2PH	U Plan Semi Apochromat phase objective 10X/0.3, WD 10		1
LUCPLFLN20XPH	Long working distance plan semiapochromat phase objective 20X, NA 0.45, WD 6.6-7.8 (c.c.0-2)		1
LUCPLFLN40XPH	Long working distance plan semiapochromat phase objective 40X, NA 0.6, WD 3.0-4.2 (c.c.0-2)		1
目镜			
WHN10X	Widefield eyepiece 10X		1
WHN10X-H	Widefield eyepiece 10X, focusable		1
对中望远镜			
U-CT30	Centering telescope		1
绿色滤光镜			
43IF550-W45	Interference green contrast filter 45mm dia.		1