



CKX41-A32PH 详细配置清单

E/C Code	Description	quantity
CKX41-A32PH 显微镜镜体	Inverted Microscope (Made in Philippines)	1
CKX41SF 三目镜筒	CKX41 frame, including 30W lampsocket AR3553, insert plate AQ1180, allen wrench GA3720, tag AD0764, 6V30WHAL (2pcs) and AB4X10SA screw Corresponded to WEEE regulation	1
U-CTR30-2 载物台延长板	Trinocular tube, FN20, fixed light pass, Bi/Photo: 50/50	1
CK2-SS 机械式载物台	Stage extension plate	1
CK40-MVR 玻璃插板	Attachable mechanical stage with right-hand vertical low drive controls (including Petri-dish holder, Terasaki- plate holder, and Slide glass holder)	1
CK40-CPG30 相差插板	Glass stage insert plate for IX50/70, CK40F and CKX41SF	1
IX2-SL 相差环板	Phase contrast slider, centrable, with phase ring for 4X, phase ring for 10X/20X and hole for phase ring for 40X (IX2-SLPH2) or filter (45mm dia.), including centering tools AQ2302 (2 pcs)	1
IX2-SLPH2 IX2-SLPHC 电源线	Phase contrast ring slit	1
UYCP 防尘罩	Phase contrast ring slit	1
COVER-012 万能平场半复消色差 相差物镜	Power cord	1
UPLFLN4XPH 平场消色差相差物镜	Dust cover Type 012 for CKX, CK (L 300mm x W 440mm x H 440mm)	1
CPLN10XPH 长工作距离消色差相 差物镜	U plan semi apochromat phase objective 4X/0.13, WD 17	1
LCACHN20XPH	Plan achromat phase objective 10X, NA 0.25, WD 10	1
LCACHN40XPH	Long working distance achromat phase objective 20X/0.4, WD 3.2	1
	Long working distance achromat phase objective 40X/0.55, WD 2.2	1

宽视野目镜		
WHB10X	Wide field eyepiece 10X, FN20	2
对中望远镜		
CT-5	Centering telescope	1
日光平衡滤色片		
45-LBD-IF	Interference light balance daylight filter, 45mm dia.	1
绿色滤色片		
43IF550-W45	Interference green contrast filter 45mm dia.	1
中性灰度滤色片		
43ND-25W45	Neutral density filter, 45mm dia.	1

---