



SZX7-1063 详细配置清单

| <b>E/C Code</b>              | <b>Description</b>  | <b>quantity</b> |
|------------------------------|---|-----------------|
| SZX7-1063<br>显微镜镜体           | Stereo microscope   | 1               |
| SZX-ZB7<br>45度倾斜双目观察筒        | 7:1 zoom microscope body  | 1               |
| SZX-BI45<br>透返两用液晶照底座        | 45 degrees inclined binocular observation tube  | 1               |
| SZ2-ILST<br>电源线              | Transmitted and reflected light LED illuminator stand<br>Corresponded to RoHS, WEEE regulations | 1               |
| UYCP<br>目镜                   | Power cord  | 1               |
| WHSZ10X-H<br>无畸变平场复消色差<br>物镜 | Eyepiece 10X with ESD capability, F.N. 22, focusable  | 2               |
| DFPLAPO1X-4                  | Distortion free plan apochromat objective 1X, W.D.81mm  | 1               |