

BX43-32P01-FLB1 详细配置清单

E/C Code	Description	quantity
BX43机架		
BX43F	Microscope frame for transmitted microscopy, including hexagonal driver ABxxxx, immersion oil 8cc Corresponded to RoHS, WEEE regulation	1
三目观察筒		
U-TR30-2	Trinocular tube	1
一般型5孔转盘		
U-5RE-2	Quintuple Revolving Nosepiece	1
右手低位载物台		
U-SVRB-4	Ceramic surface mechanical stage with right-hand low drivecontrol (long type)	1
30瓦卤素灯适配器		
U-LS30ADP	Adapter for attaching LS30 to BX43/43C frame, with connection cable Corresponded to RoHS, WEEE regulations	1
左手双标厚型标本夹		
U-HLDT-4	Specimen holder for BX stage, left hand, for two specimens, thick type	1
摇入摇出型聚光镜		
U-SC3	Swing-out condenser N.A. 0.9-0.16	1
30瓦卤素灯电源		
TL-4	Power supply for 30W halogen lamp with output cable (0.64m) (Made in Philippines) Corresponded to WEEE regulation	1
30瓦卤素灯灯座		
U-LS30-5	Lamp socket	1
6伏30瓦卤素灯		
6V30WHAL	30W Halogen bulb (Philips No. 5761)	2
日光平衡滤色片		
45-LBD-IF	Interference light balance daylight filter, 45mm dia.	1
电源线		
UYCP	Power cord	2
防尘罩018		
COVER-018	Dust cover Type 018 for BX2 (L 510mm x W 320mm x H 530 mm)	1
平场消色差物镜		
PLN4X	Plan achromat objective 4X/0.1, WD 18.5	1

PLN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLN20X	Plan achromat objective 20X/0.4, WD 1.2 (spring)	1
PLN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
PLN100XO	Plan achromat objective 100X/1.25, WD 0.15 (spring, oil)	1
10倍宽视场目镜		
WHN10X	Widefield eyepiece 10X	1
10倍宽视场可调目镜		
WHN10X-H	Widefield eyepiece 10X	1
8毫升镜油		
OIL-8CC	Immersion oil 8cc	1
荧光装置		
中性灰度滤色片		
BX3-25ND25	Neutral density filter for BX3 incident light illuminator	1
BX3-25ND6	Neutral density filter for BX3 incident light illuminator	1
100汞灯灯室		
U-LH100HG	Lamp house for 100w mercury	1
汞灯电源		
U-RFL-T	Power supply unit for Mercury lamp meeting the RoHS regulation	1
电源线		
UYCP	Power cord	1
荧光激发块		
U-FUW	Mirror Unit for fluorescence microscopy, wide band UV excitation, exciter filter BP330-385, dichroic beamsplitter DM400, barrier filter BA420	1
U-FBW	Mirror unit for fluorescence microscopy, wide band blue excitation,exciter BP450-480,dichroic beamsplitter DM500,barrier filter BA515	1
U-FGW	Mirror unit for fluorescence microscopy, wide band interference green excitation, with exciter filter BP530-550, dichroic beamsplitter DM570 and barrier filter BA575-UF, including plastic case ARxxxx	1
BX3通用荧光反射器		
BX3-URF	Universal reflected light illuminator with field stop and aperture stop, including 8-position mirror unit cassette, UV protective shield ADxxxx, shield plate AB1550,Allen wrench AD6162 and indication plate ADxxxx	1
荧光对中测试标本		
U-CST	Centering specimen with target for fluorescence microscopy	1
100瓦汞灯灯泡		
USH103OL	100W mercury burner for U-LH100HGAP0,U-LH100HG	

从普通，到声名显赫，普赫 (Pooher) 都未离开您 电话：021-61678408 www.pooher.com

30毫升无荧光观察用镜油

IMMOIL-F30CC

Immersion oil 30cc. Low level auto-fluorescence.

1
