

### BX53T-32P02-FLB1详细配置清单

E/C Code	Description	Q
BX53机架	Microscope frame for transmitted and reflected microscopy with built-in filter holders for four filters (LBD, ND6, ND25 and GF), including hexagonal driver ABxxxx and immersion oil 8cc Corresponded to WEEE and RoHS regulations	1
BX53F		
三目观察筒	Trinocular tube	1
U-TR30-2		
D型6孔物镜转盘	Sextuple revolving nosepiece with a slot for analyzer or DICslider	1
U-D6RE		
右手低位载物台	Ceramic surface mechanical stage with right-hand low drivecontrol (long type)	1
U-SVRB-4		
左手双标厚型标本夹	Specimen holder for BX stage, left hand, for two specimens, thick type	1
U-HLDT-4		
摇入摇出型聚光镜	Swing-out condenser N.A. 0.9-0.16	1
U-SC3		
100W 卤素灯灯室	Lamp house for 100W halogen	1
U-LH100-3		
12伏100瓦长寿命卤素灯泡	Halogen bulb (Philips No. 7724)	2
JC12V100WHAL-L		
电源线	Power cord	1
UYCP		
BX3镜臂	Arm for BX53F and BX63F frame	1
BX3-ARM		
防尘罩018	Dust cover Type 018 for BX2 ( L 510mm x W 320mm x H 530 mm)	1
COVER-018		
平场消色差物镜	Plan achromat objective 4X/0.1, WD 18.5	1
PLN4X		
PLN10X	Plan achromat objective 10X/0.25, WD 10.6	1
PLN20X	Plan achromat objective 20X/0.4, WD 1.2 (spring)	1
PLN40X	Plan achromat objective 40X/0.65, WD 0.6 (spring)	1
PLN100XO	Plan achromat objective 100X/1.25, WD 0.15 (spring, oil)	1
10倍宽视场目镜		

WHN10X 10倍宽视场可调目镜	Widefield eyepiece 10X	1
WHN10X-H 中性灰度滤色片	Widefield eyepiece 10X, focusable	1
BX3-25ND25	Neutral density filter for BX3 incident light illuminator	1
BX3-25ND6	Neutral density filter for BX3 incident light illuminator	1
100汞灯灯室		
U-LH100HG 汞灯电源	Lamp house for 100w mercury	1
U-RFL-T 电源线	Power supply unit for Mercury lamp meeting the RoHS regulation	1
UYCP 荧光激发块	Power cord	1
U-FUW	Mirror Unit for fluorescence microscopy, wide band UV excitation, exciter filter BP330-385, dichroic beamsplitter DM400, barrier filter BA420	1
U-FBW	Mirror unit for fluorescence microscopy, wide band blue excitation, exciter BP450-480, dichroic beamsplitter DM500, barrier filter BA515	1
U-FGW BX3通用荧光反射器	Mirror unit for fluorescence microscopy, wide band interference green excitation, with exciter filter BP530-550, dichroic beamsplitter DM570 and barrier filter BA575-UF, including plastic case ARxxxx	1
BX3-URF	Universal reflected light illuminator with field stop and aperture stop, including 8-position mirror unit cassette, UV protective shield ADxxxx, shield plate AB1550, Allen wrench AD6162 and indication plate ADxxxx	1
荧光对中测试标本		
U-CST 100瓦汞灯灯泡	Centering specimen with target for fluorescence microscopy	1
USH103OL 30毫升无荧光观察用 镜油	100W mercury burner for U-LH100HGAP0,U-LH100HG	
IMMOIL-F30CC	Immersion oil 30cc. Low level auto-fluorescence.	1