



CX22LEDRFS1 详细配置清单

Units	Description	Qty
教学级生物显微镜		
CX22LEDRFS1	Biological Microscope	1
CX22 显微镜镜体 (配置 4X/10X/40X/100X 物镜)	Microscope frame for transmitted microscopy with binocular tube, a pair of eyepiece 10X (F.N.20) , quadruple revolving nosepiece, mechanical stage, abbe condenser, plan objective (4X, 10X, 40X, 100X) , 0.5W LED Internal	1
CX22LEDRFS1		1
电源线		
UYCP	Power cord	1
防尘罩		
COVER-015	Dust cover Type 015 for CX22/CX21/31/41 (L 450mm x H 460mm)	1