

IX73-F22FL/PH 详细配置清单

Units	Description	Qty
显微镜镜体		
IX73P1F	IX73P1F manual frame with 1port for fluorescence illuminator/ magnification changer/ right port, etc. 3 position light path prism (100:0/50:50/0:100) integrated. Motorized/coded components can be attached with IX3-CBM/U-CB5S. Light intensity operation and	1
编码型 6 孔物镜 转盘		
IX3-D6RES	IX3-D6RES Coded 6 position nosepiece for IX73P2F/P1F. DIC slider slot is available. U-CB5S is necessary to utilize coded function. Combination with transmitted LED (U-LEDPS) enable to use LIGHT MANAGER function. U-NPCBL100 is required for U- CB5S/U-LEDPS	1
双目镜筒		
U-BI90	Binocular tube	1
目镜		
WHN10X-H	Widefield eyepiece 10X, focusable Made in China	1
WHN10X	Widefield eyepiece 10X Made in China	1
对中望远镜		
U-CT30-2	Centering telescope Corresponded to RoHS regulation	1
卤素灯照明支架		
IX3-ILL	IX3-ILL transmitted pillar for IX83P2ZF/P1ZF, IX73P2F/P1F, IX53F. Front operation knob for condenser focus, focus stopper for reproduce position. Halogen and LED light source can be mounted. (Made in China).	1
卤素灯电源		
TH4(200V)	Power supply	1
卤素灯灯箱		
U-LH100L-3	Lamp house for 100W halogen including long connection cord Corresponded to RoHS, WEEE regulations	1
延长电源线		
U-RMT	Extension cable for halogen lamphousing	1
卤素灯灯泡		
JC12V100WHAL-L	Halogen bulb (Philips No. 7724)for AX70/80,BX50/51/52/60/61/62IX70 /71/81,AHBS/3,AHBT/3,BHS/U IMT2 (100W) (Made in Germany)	2
滤波滤色片		
45FR	Frost filer for IX3-ILL. Diameter is 45mm.	1
日光平衡滤色片		

45-LBD-IF 减光滤色片	Interference light balance daylight filter, 45mm dia.	1
45-ND6	Neutral density filter, 45mm dia.	1
45-ND25	Neutral density filter, 45mm dia.	1
连接电源线		
UYCP-33 机械式载物台	Power cord	1
IX3-SVR	IX3-SVR mechanical right hand stage with fixed handle. Enough travelling range applied for slide glass, 35mm dish as well as multiwell plate. circle stage inserts included. Sample holder for slide glass, 35mm dish, multiwell plate as option. Stage stopper	1
玻璃载物台中心板		
IX2-GCP	Glass insert plate with 25mm dia. Opening, including rubber bands for objective lens (5 pcs; Red (AD9322), Yellow (AD9323), Green (AD9324), Blue (AD9325), White (AD9326), Light Blue (AD9327))	1
长工作距离万能聚 光镜		
IX2-LWUCD	Long working distance universal condenser (N.A.0.55/W.D.27mm) with 5 position turret (three for 30mm dia. And two for 38mm dia.) and aperture stop, including centering tools AQ4512 (2 pcs) and indication plate AC1338	1
相差环板		
IX-PHL	Phase contrast ring slit (30mm dia.)	1
IX-PH1	Phase contrast ring slit (30mm dia.)	1
IX-PH2	Phase contrast ring slit (30mm dia.)	1
绿色滤光片		
43IF550-W45 万能平场半复消色 差相差物镜	Interference green contrast filter 45mm dia.	1
UPLFLN4XPH	U plan semi apochromat phase objective 4X/0.13, WD 17	1
UPLFLN10X2PH	U Plan Semi Apochromat phase objective 10X/0.3, WD 10	1
长工作距离万能平 场半复消色差相差 物镜		
LUCPLFLN20XPH	Long working distance plan semiapochromat phase objective 20X, NA 0.45, WD 6.6-7.8 (c.c.0-2)	1
LUCPLFLN40XPH	Long working distance plan semiapochromat phase objective 40X, NA 0.6, WD 3.0-4.2 (c.c.0-2)	1
编码型落射荧光激 发块转盘		
IX3-RFACS	IX3-RFACS Coded fluorescence mirror turret with 8 Position. Manual shutter is available. U-CB5S/ U-COCBL100 are	1

	required to utilize coded function to read out position.	
落射荧光照明器 IX3-RFA	IX3-RFA straight type fluorescence illuminator. Field stop is equipped. 1 slot for neutral density filter. Illuminator is designed and optimized for upper deck for IX83P2ZF/IX73P1F, and IX83P1F/IX73P1F. (Made in China).	1
落射光路减光滤色片 25ND6	Neutral density filter, 6 percent transmission, 25mm dia.	1
25ND25	Neutral density filter, 25 percent transmission, 25mm dia.	1
汞灯灯箱 U-LH100HG	Lamp house for 100W mercury burner with collector lens and connecting cable	1
汞灯灯泡 USH-103OL	100W mercury burner for U-LH100HG APO, U-LH100HG	2
汞灯电源 U-RFL-T	Power supply unit for Mercury lamp meeting the RoHS regulation	1
荧光显微镜对中标本 U-CST	Centering specimen with target for fluorescence microscopy	1
汞灯对中调节手柄 U-CLA	Collector lens adjustment knob for U-ULH	1
连接电源线 UYCP-33	-	1
荧光激发块 U-FUW	Mirror Unit for fluorescence microscopy, wide band UV excitation, exciter filter BP340-390, dichroic beamsplitter DM410, barrier filter BA420-IF, including plastic case AR1924	1
U-FBW	Mirror Unit for fluorescence microscopy, wide band Blue excitation, exciter filter BP460-495, dichroic beamsplitter DM505, barrier filter BA510-IF, including plastic case AR1924	1
U-FGW	Mirror unit for fluorescence microscopy, wide band interference green excitation, with exciter filter BP530-550, dichroic beamsplitter DM570 and barrier filter BA575-IF, including plastic case AR1924	1
防尘罩 COVER-020	Dust cover Type 020 for IX2, MVX10 (L 620mm x W 250mm x H 660 mm)	1